

## BATH/SHOWER STOOLS



42900

• **42900 SHOWER STOOL - fixed**  
Shower stool with aluminium frame and plastic seat 50x30 cm. Height adjustable 35-45 cm. Non slip end. Load: 100 kg




42901

• **42901 SHOWER STOOL - foldable**  
Shower stool with aluminium frame and plastic seat 50x30 cm. Height adjustable 44-47.5 cm. Non slip end. Load: 100 kg



42903

• **42903 SHOWER STOOL**  
Shower stool with aluminium frame and plastic seat Ø 32.5 cm. Height adjustable 38-56 cm. Non slip end. Load: 100 kg.



Adjustable seat 54.4 - 71.6 cm with 8 height adjustments, by 2.5 cm (1") increment

27733

• **27733 ADJUSTABLE SWIVEL STOOL WITH HANDLE**  
- swivel seat with 90 degree position lock at every right angle location  
- non-slip central seat pad and soft under-covers  
- load 135 kg (300 lbs)

• **42905 SHOWER STOOL WITH ARMRESTS**

Shower stool made of silver aluminium frame with PU seat, foam armrest and adjustable foot tube. Weight: 4.0 kg Load: 136 kg Size: 44.8-47.2x50-52.6x h 62.2-72 cm Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.



42905

## BATH STEPS

• **42920 BATH STEP - 1 step**  
• **42921 BATH STEP - 2 steps**  
Stainless steel frame and plastic steps. Load: 100 kg



42920

Size: 27x41xh 14 cm



42921

Size: 51x41xh14-29 cm

## SAFETY SUCTION GRAB BARS



28974




28972

28973

**VISIBLE SAFETY**  
GREEN means well attached  
RED is a warning signal. Need re-mount

• **28972 SAFETY GRAB BAR - 443 mm**  
• **28973 SAFETY GRAB BAR - 545 mm**  
• **28974 SAFETY DOUBLE GRAB BAR - 929 mm**  
Innovative line of GRAB BAR, can be quickly installed everywhere without need of drilling and screwing. Patented suction cups (Ø 93 mm). Holds up to 60 kg (28972/3) or 85 kg (28974) vertical and 45 or 60 kg (28974) lateral weight (tested by SGS). Ideal solution to support elder people during daily bathroom use. Top quality product made of heavy duty ABS, with patented "green/red" suction indicator for better safety. Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic.

## PLASTIC RIBBED GRAB BAR



28977

28976

MADE IN ITALY

• **28975 PLASTIC RIBBED GRAB BAR - 30 cm**  
• **28976 PLASTIC RIBBED GRAB BAR - 45 cm**  
• **28977 PLASTIC RIBBED GRAB BAR - 60 cm**  
Ribbed plastic grab bar in Medical PVC, complete with ABS mounting plates and screw covers. Wall screws not included. Electrically isolated. The ribbed bars offer excellent anti-slip grid and are warm to the touch. Can be mounted either horizontally or vertically. Tube Ø 3.5 cm. Max user weight 120 kg (265 lb). Manual: GB, IT.

## U-SHAPE FOLDING BAR



28979

MADE IN ITALY

• **28979 U-SHAPE POWDER COATED STEEL FOLDING BAR - 70 cm**  
White powder coated wall mounted lift up toilet support. Designed to help getting up or transferring to and from the toilet. It can be folded upwards for easy accessibility. Max user weight 125 kg (275 lb). Wall screws not included. Manual: GB, FR, IT, ES, DE, GR.

